



Using a controlled vacuum

V-ITEX

VACUUM IN TUBE EXTRACTION

Automated vacuum extraction for GC-MS analysis, simplified

V-ITEX

Vacuum In Tube Extraction


GERSTEL V-ITEX is an innovative, fully automated, vacuum-based extraction method developed by Agroscope (Switzerland) for determining volatile organic compounds (VOCs) in complex matrices using GC-MS. Using a controlled vacuum, V-ITEX enables gentle extraction and enrichment even from temperature-sensitive and complex matrices – without heating, with virtually no solvent consumption, and with maximum reproducibility.

After enrichment, the analytes are directly transferred to the GC-MS system, providing extremely low detection limits with high precision and accuracy.

Challenges of VOC extraction

Many laboratories face the same problems when extracting VOCs:

- **Time-consuming and labor-intensive methods**
such as SAFE or Likens-Nickerson demand valuable resources
- **High solvent consumption**
incurs costs and jeopardizes occupational safety in the laboratory
- **Non-volatile components**
cannot be analyzed using conventional headspace GC-MS methods
- **Artifact formation and falsified results**
due to the heating of temperature-sensitive samples



Features


- Gentle VOC extraction of thermally labile analytes using vacuum and no heat
- Automatic enrichment and water management with the GERSTEL MultiPurposeSampler (MPS) and the V-ITEX module
- Direct desorption of the trap utilizes GC carrier gas stream

Areas of application

- Flavor analysis in dairy products, beverages, fermented foods
- Nutritional volatolomics VOC profiles in food and biological fluids
- Temperature-sensitive foods such as fresh food, milk, soy milk, and oat milk
- Research in the field of volatolomics and VOC metabolomics

Highlights

- Efficient drying enables total evaporation even of aqueous samples for improved recovery of polar compounds
- Simple and robust automation for unattended, worry-free operation
- GERSTEL Agitator Stirrer (AGI^{stir}) enables parallel incubation of multiple samples for higher throughput
- Vacuum-controlled pump (Vacuubrand) for reproducible results
- Seamless integration into almost any GC or GC-MS through GERSTEL ePneumatics
- Agilent MassHunter Acquisition Software Integration



Automated vacuum extraction for GC-MS analysis, simplified



GERSTEL, Inc.

Japan
GERSTEL K. K.
 ☎ + 81 3 5731 5321
 ✉ info@gerstel.co.jp

📍 **Eberhard Gerstel Platz 1**
 45473 Mülheim an der Ruhr
 Germany

📍 **701 Digital Drive, Suite J**
 Linthicum, MD 21090
 United States

Singapore
GERSTEL LLP
 ☎ +65 6779 0933
 ✉ sea@gerstel.com

☎ +49 208 76503-0
 ✉ info@gerstel.com
 🌐 **www.gerstel.com**

☎ 410.247.5885
 ✉ sales@gerstelus.com
 🌐 **www.gerstelus.com**

China
GERSTEL (Shanghai) Co. Ltd
 ☎ +86 (0) 21 507 19 39
 ✉ china@gerstel.com

Velaris reserves the right to change the specifications and the appearance of the equipment without further notification.

Version 1.1